

WIR RX15-2 Receiver - 2.3 MHz or 2.8 MHz

INSTRUCTIONS FOR USE AND CARE

Safety Information

HEARING SAFETY:

CAUTION!

This product is designed to amplify sounds to a high volume level which could potentially hearing damage if used improperly. To protect your hearing and the hearing of others:

1. Make sure the volume is turned down before putting on the earphone or headphone before adjusting the volume to a comfortable level.
2. Set the volume level at the minimum setting that you need to hear.
3. Do not allow children or other unauthorized persons to have access to this product.

BATTERY SAFETY & DISPOSAL

CAUTION!


Do not attempt to recharge disposable batteries, which may explode, release dangerous chemicals, cause burns, or other serious harm to the user or product.

PACEMAKER SAFETY

CAUTION!

1. Before using this product with a pacemaker or other medical device, consult your physician or the manufacturer of your pacemaker or other medical device.
2. If you have a pacemaker or other medical device, make sure that you are using this product in accordance with safety guidelines established by your physician or the pacemaker manufacturer.

Operating Instructions



A127

Battery Instructions

Two batteries (AAA Alkaline) must be placed into the battery compartment of the headphone, and ensure that the positive and negative terminal on the batteries are inserted correctly.

Operation Instructions

1. Slide the Power and channel switch of the headphone to the position A / B to power on the headphone. The Power and channel indicator turns red when the Power and channel switch slides to A and turns green when the Power and channel switch slides to B. You can enjoy different varieties of audio with different channel selections.
2. Adjust the volume of the headphone with the Volume control knob to desired level.
3. You can silence the audio to hear outside environment by holding the Monitor button continuously and releasing the Monitor button to continue listening to audio on selected channel.
4. Slide the Power switch of the headphone to position OFF to power off the headphone when it is not in use.
5. Connect audio source and Audio input jack of the headphone with audio connect line, the headphone is also used as regular wired stereo headphone without battery.

Special Features

1. The WIR RX15-2 receiver headphone has an auto mute function: When the IR transmitter turns off or the WIR RX15-2 does not receive the IR signal, the WIR RX15-2 automatically turns off its power amplifier and mutes.
2. The WIR RX15-2 receiver headphone has an auto power off function: After 3-5 minutes if the WIR RX15-2 does not receive the IR signal, the WIR RX15-2 will automatically turn off its power and the power indicator will turn off. To reactivate when the signal is present, the WIR RX15-2 power switch should be turned off and then back on.
3. The WIR RX15-2 headphones are equipped with a moving coil type speaker and are compatible with T-coil equipped hearing aids.
4. Audio input jack. Use optional cord to plug into audio source.

Specifications

WIR RX15-2 Receiver

System	Infrared light (850nm)
Reception	360 degrees
Carrier Frequency	Channel A: 2.3MHz
	Channel B: 2.8MHz
Modulation	FM
Mode	MONO
S/N	> = 60dB
Distortion	< 1.5% (1KHz sine wave)
Max Average SPL	118dB
Frequency Response	30-20,000Hz
Power	(2) AAA Alkaline batteries
Power Indicator	Red = 2.3MHz; Green = 2.8MHz
Battery Life	Up to 80 hours
Current Draw	Normal 12mA
Approvals	CE, RoHS, WEEE

Recycling Instructions

Help Williams Sound protect the environment! Please take the time to dispose of your equipment properly.



Product Recycling:

Please do NOT dispose of your Williams Sound equipment in the household trash. Please take the equipment to a electronics recycling center; OR, return the product to the factory for proper disposal.



Battery Recycling:

Please do NOT dispose of used batteries in the household trash. Please take the batteries to a retail or community collection point for recycling.

Approvals

This product meets approval standards established for:

CE RoHS WEEE

Limited Warranty

Williams Sound products are engineered, designed and manufactured under carefully controlled conditions to provide you with many years of reliable service. Williams Sound warrants the Infrared Receiver against defects in materials and workmanship for TWO (2) years. During the first two years from the purchase date, we will promptly repair or replace the Infrared Receiver.

Earphones, headphones, batteries, cables, carry cases, and all other accessory products carry a 90-day warranty.

WILLIAMS SOUND HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. WILLIAMS SOUND, THEREFORE, DISCLAIMS ALL WARRANTIES NOT SET FORTH ABOVE, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE INFRARED RECEIVER, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. WILLIAMS SOUND SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. THE SOLE REMEDY FOR ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT IS REPLACEMENT OF THE PRODUCT. NO PERSON HAS ANY AUTHORITY TO BIND WILLIAMS SOUND TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE INFRARED RECEIVER. UNAUTHORIZED REPAIRS OR MODIFICATIONS WILL VOID THE WARRANTY.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

Customer Service

If you experience difficulty with your system, call Toll-Free for Customer Assistance:

1-800-843-3544

If it is necessary to return the system for service, your Customer Service Representative will give you a Return Authorization Number (RA) and shipping instructions.

Pack the system carefully and send it to:

Williams Sound
Attn: Repair Dept.
10300 Valley View Rd
Eden Prairie, MN 55344



10300 Valley View Rd • Eden Prairie, MN 55344
800-328-6190 / 952-943-2252 • FAX: 952-943-2174
www.williamssound.com

Printed in the USA
©2012 Williams Sound, LLC
All Rights Reserved

MAN162B